

Vortex / Cyclomixers

Vortex Mixers

- Designed for mixing liquids in Schools, Laboratories & Factories
- Touch/ Continuous Operation mode Selection through bi-directional Switch
- Speed Regulation through knob provided on the control panel
- Interchangeable mixing heads for use with variety of tubes



CM - 101 Plus


Technical Data	
Models	CM - 101 PLUS
Power [W]	60
Shaking Movement	Orbital
Orbital Diameter [mm]	4
Motor Type	Shaded-Pole Motor
Motor Rating Input / Output [W]	58 / 10
Permissible ON time	100% power 30mins
Speed range [rpm]	0-2500
Run Type	Continuous / touch operation
Dimensions (mm)	127 x 130 x160
Weight [kg]	3.5
Protection Class acc. To DIN 60529	IP 21

Supply 220/240 Volts, 50 Hz


Interchangeable Mixing Heads	
Models	Description
CM-01	Universal top plate 100 mm-plate attachment
CM-01/01	Tube adapter for 18 holes test tubes ø 10mm (use with CM-01)
CM-01/02	Platform pad for <ø99mm tubes and small vessels (use with CM-01)
CM-01/03	Tube adapter for 8 holes test tubes ø 20mm (use with CM-01)



CM - 01



CM - 01/ 03



CM - 01/ 01